

SCANNER	HOW TO SEND
TRIOS	<ol style="list-style-type: none"> Go to www.us.3shapecommunicate.com and enter your Username and Password. Click on “Connections”, then “Add Connection”. Enter: melissa@midwestdentalarts.com We will be notified to approve your connection and once this is done, you will be able to send us cases. <p><i>If any issues arise, please contact Henry Schein at 800-537-6070</i></p>
CARESTREAM (CS3500)	<ol style="list-style-type: none"> After acquiring scan data, you have two options to send us files, CSConnect and Send to: <ol style="list-style-type: none"> CS Connect: Login w/ email & password, choose our name under “pick lab” or input our lab’s email address if this is the first case: cad@midwestdentalarts.com Send to: gives you the option of saving the file to the desktop where you could then attach the files in an email to cad@midwestdentalarts.com. <p><i>If any issues arise, please contact Carestream tech support at 866-724-6317 opt 6 and 1 or 2</i></p>
CEREC CONNECT	<ol style="list-style-type: none"> Go to cerecconnect.com, login, add MidwestDentalArts to your laboratory list. Search by name (Midwest Dental Arts) and location Swisher/Cedar Rapids, IA Open CEREC CONNECT SW4 Icon To start a case: enter patient information, tooth number, restoration type, and shade Scan patient and draw margins Review the scan quality, interocclusal clearance, and margins. Select “NEXT”- this button takes you to the CEREC CONNECT PORTAL Choose Midwest Dental Arts as laboratory and Send case. <p><i>You will receive an email notification that the file was sent. If any issues arise, please call CEREC SUPPORT at 800-475-5036</i></p>
iTERO	<ol style="list-style-type: none"> Create New Case Fill out Treatment Information (patient info, scanning order, restoration type, and shade). Choose Midwest Dental Arts as lab of choice. To set up Midwest Dental Arts in your system, call iTero customer support if we are not in your system. Our acct # is: 2011 Scan Patient Review the scan quality, interocclusal clearance, and margins. Click on the “Send” icon in the main toolbar to send the case to the Cadent Center. To confirm that the file has been sent (or is queued for sending) open the Case Manager. <p><i>If any issues arise, please contact iTero Customer Service at 800-577-8767</i></p>
LAVA COS	<ol style="list-style-type: none"> Add Midwest Dental Arts to your lab list by calling 3M Customer Support at 800-634-2249; ext 3, ext 1 for digital tech support. After the lab is set up, it must also be activated. Enter new patient information Scan patient Review the scan quality, interocclusal clearance, and margins. Select Midwest Dental Arts as lab and select material for restoration Sign “Prescription”. <p><i>Contact Customer Support, should issues arise at 800-634-2249, ext 3, ext 1</i></p>
TRUDEF	<p>Add Midwest Dental Arts to your lab list by calling 3M Customer Support at 800-634-2249; ext 3, ext 1 for digital tech support. After the lab is set up, it must also be activated.</p> <ol style="list-style-type: none"> Enter new patient information Scan patient Review the scan quality, interocclusal clearance, and margins. Select Midwest Dental Arts as lab, <u>No Models</u>, and select material for restoration <p><i>Contact Customer Support, should issues arise at 800-634-2249, ext 3, ext 1</i></p>
E4D	<ol style="list-style-type: none"> You must be registered through DDX in order to send cases to us. Contact Henry Schein for technical assistance to add Midwest Dental Arts, Inc to your list of laboratories: 800-537-6070.